



### STANDARD SIZES AVAILABLE

S.No.	Sizes Available (diaXHeight)	Cap. in Ltr.
1.	25 x 45 cm	22
2.	30 x 50 cm	35
3.	35 x 55 cm	52
4.	45 x 60 cm	95
5.	55 x 75 cm	175
6.	55 x 90 cm	200

## VERTICAL AUTOCLAVE

- Reliable Convenient Method Of Steam Sterilization.
- Inner Chamber Made Of Thick S. S.
- Outer Chamber Made Of M.S / S. S. Immersion Type Heater Should Be Provided. Lid Is Made Up Of S. S. and Die Pressed.
- It Is Equipped With Pressure Gauge, Steam Release Valve And Safety Valves.
- This Instrument Supplied With S. S. Basket.
- Pressure Controls By Spring Valves.
- It Is Used For Sterilization Under Working Steam Pressure Upto 15 Psi.
- Fitted With Silicon Rubber Gasket Jointless, To Work On 230 Volts A./C. Only.
- Foot Lifting Arrangement.

#### Optional accessories:

- A) Foot Lifting Arrangement.
- B) Digital Temperature Indicator Cum Controller with Timer.
- C) Water Level Indicator.
- D) Pressure cutoff Switch.



## SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : [www.lutrontestinginstruments.com](http://www.lutrontestinginstruments.com)