



STANDARD SIZES AVAILABLE		
S.No.	Sizes Available (diaXHeight)	Cap. in Ltr.
1.	25 x 45 cm	22
2.	30 x 50 cm	35
3.	35 x 55 cm	52
4.	45 x 60 cm	95
5.	55 x 75 cm	175
6.	55 x 90 cm	200

VERTICAL AUTOCLAVE FULLY AUTOMATIC (BTI-02-A)

- Tripal Body Construction Separate Chamber For Steam & Water Covered Externally By Outer Shell.
- Lid Made Of Thick Stainless Steel Pressure Up To 22 Psi Factory Set At 15 Psi (Temp 121°C)
- Complete Unit With Automatic Air Purging Digital Temp. Controller & Indication With 0.1 C Accuracy/ Resolution With Electronic Digital Timer Range 0 To 99 Min
- Autoclaving Time Cycle Automatically Starts When Set Pressure Is Achieved.
- Automatic Steam Release Once Autoclaving Cycle Complete.
- Complete With Pressure Gauge, Safety Valve Stem Release Valve, And Extra Safety Valve.
- Low Water Level Cut-off: Low Water Level Cut-off Assembly Connected To Heater Which Automatically Cut-off Power To Heater If Water Level Is Not Sufficient Thus Eliminating Heater Burn Out.
- Foot Lifting Arrangement For Lid.
- Provided With S. S. 304 Grade Perforated Basket.

Manufactured by :

BIO - TECHNICS INDIA™

AN ISO 9001 : 2008