



**IMPORTERS / EXPORTERS & SUPPLIERS :-**

# **SUNSHINE INSTRUMENTS**

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com

## **BULK DENSITY APPARATUS**



### **BULK DENSITY APPARATUS (DIGITAL)**

Microprocessor based 4 digit preset counter 30 strokes per minutes, strokes height 50 mm(Supplied without cylinder)

### **BULK DENSITY APPARATUS (ANALOGUE)**

With 15 Min Mechanical Timer 30 strokes per minutes, stroke height 50 mm (supplied without cylinder)

# **SUNSHINE INSTRUMENTS**