



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com

BALL MIL 1 KG RBMT-1 /R BMT-2



A drum is made of welded steel having an inside diameter of 300 mm x 300 mm long. Supported on heavy duty ball bearings. An air and watertight cover is provided to close the 100 mm. Wide opening along with entire length of the drum. A reduction gear driven by a motor rotates the mill at 28 - 30 rpm. Supplied complete with twelve 19 mm. Diameter hardened steel balls and an automatic revolution counter. Operation on 1 H.P. motor, 220 / 230 V, 50 Hz. Single Phase Main AC.

Available in difference sizes.

- 1.) Ball Mill 1 Kg (With Aluminium Jar)
- 2.) Ball Mill 1 Kg (With S.S. Jar)
- 3.) Ball Mill 5 Kg. (in Mild Steel Drum)
- 4.) Ball Mill 5 Kg. (In S.S. Drum).