



**IMPORTERS / EXPORTERS & SUPPLIERS :-**

# SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com

## BUNSEN BURNER



### BUNSEN BURNER: ELECTRIC

- ☞ This apparatus consists of heating element.
- ☞ It is in conically shaped refractory crucible.
- ☞ It is enclosed in an insulated S.S. casing.
- ☞ Casing is mounted on top of S.S. Cylindrical housing.
- ☞ Maximum working temp. 700°C. approximately
- ☞ Temp. controlled by energy regulator.
- ☞ To work at 230 Volts AC Only.

**BUNSEN BURNER:** With temperature controller, Maximum Temp.700°C

Supplied with body of base made of M.S Powder coated & top body made of S.S with energy regulator to control temperature. Maximum surface temperature upto 700° C

SPARE CRUCIBLE FOR BUNSEN BURNER

# SUNSHINE INSTRUMENTS