



**IMPORTERS / EXPORTERS & SUPPLIERS :-**

## **SUNSHINE INSTRUMENTS**

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com

### **FRIABILITY TEST APPARTUS (Non Digital)**



### **FRIABILITY TEST APPARTUS (Non Digital)**

Mechanical timer for 5 Minutes.25 revolutions per 1 Min .universal acrylic drum for right or left side easy loading & unloading of drum. Provide with drum

### **FRIABILITY TEST APPARTUS (DIGITAL )**

Digital Timer preset counter Upto 99.9 minutes,25 revolution per 1 min universal acrylic drum for right or left side easy loading & unloading of the drum provide with drum.

# **SUNSHINE INSTRUMENTS**