



**IMPORTERS / EXPORTERS & SUPPLIERS :-**

# SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : [instrumentsindia2012@gmail.com](mailto:instrumentsindia2012@gmail.com)

Website : [www.lutrontestinginstruments.com](http://www.lutrontestinginstruments.com)

## **FILTER HOLDER: 47mm S.S Filter Holder**

### **Specification:**

Material: Stainless steel 316 Quality

Filter support: stainless steel screen

Sealing: silicon O -ring

Funnel Cap: 250 MI

Filter diameter required: Membrane filter, 47 mm pre filter, and 47mm

Filtration area: 12.5 cm<sup>2</sup>

Maximum operating pressure: Vacuum Only.

Sterilization: Autoclaving (max.121 °C)

**ALSO AVAILABE in 13 mm / 25 mm**



# SUNSHINE INSTRUMENTS

 [instrumentsindia2012@gmail.com](mailto:instrumentsindia2012@gmail.com)

 +91 93600 23357