



TX223L, TX323L  
TX423L, TXC

## TXC Series

### TXC series

Model	TXC323L	TXC623L
Capacity	320ct	620ct
Minimum Display	0.001ct	
Repeatability	0.001ct	
Linearity	±0.002ct	
Operating Ambient Temperature	5-40°C	
Temperature Coefficient Sensitivity (ppm/C)	±3	
Pan Size (mm) approx.	ø80	
Main Body Dimensions (mm) approx.	206 ( W ) × 291 ( D ) × 241 ( H )	
Weight (kg) approx.	3.8kg	
Power requirement	AC adaptor ( DC12V, 1A )	

