



IMPORTERS / EXPORTERS & SUPPLIERS :-

# SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com



## THERMOMETER GLASS ZEAL 300mm

Zeal - 10c to 50°C mercury filled Size 300mm Accuracy 0.5°C

Zeal - 10c to 110°C (RED ALCHOL FILLED) Size 300mm Accuracy 1°C

Zeal - 10c to 150°C (M) Size 300mm Accuracy 1°C

Zeal -10 to +50 or 110 or 150 Alcohol filled Size 300mm Accuracy 1°C (Anyone)

Zeal - 0c to 250°C (M) Size 300mm Accuracy 1°C

Zeal - 0c to 360°C (M) Size 300mm Accuracy 2°C

Zeal - 0c to 400°C (M) Size 300mm Accuracy 2°C

Zeal 0c - 105°C (1/5th) (M) Size 375mm Accuracy 0.2°C

Zeal 0c - 50°C (1/10th) (M) Size 375mm Accuracy 0.1°C

zeal pocket -20 to +20F Thermometer (A) Size 150mm Accuracy 2F

zeal - 0 to 150F / 240F Mercury Size 300mm Accuracy 1F

Zeal -400F Size 300mm Accuracy 2F

Spare Brass cover for 300mm Stem Glass Thermometer, Indian

Certification for Thermometer extra. (Traceability to National standard) Issued by private party

Zeal Curo Meter : Tobacco Barm Thermometer

Zeal -50°C + 50°C Alchol Toulene Filled - Red

