




DIGITAL FRIABILITY TEST APPARATUS

This apparatus is designed for determination of abrasion and impact hardness of the tablets.

SALIENT Features

- ❖ Programmable Rotation Time with digital display
- ❖ Motorized Addition/Removal of tablets without opening the drum.
- ❖ Auto/Manual mode of operation.
- ❖ Special Low RPM motor means less noise and constant speed even at very Low/High Voltage (160-260).
- ❖ Auto Stop with audio alarm.



TECHNICAL DETAILS

Description	DVSFA-263	DVSMP-264
Rotation Speed	25+-1rpm	25+-1rpm
Display	3 digit LED	3digit LED
No of Drums	One	Two
Electric Motor	Low rpm A.C.Motor	Low rpm A.C.motor
Audio Alarm	Yes	Yes
Set Time	Programmable	Programmable

Due to Continuous development, look and specifications may change.

This is ERBA MAKE