



SUNSHINE INSTRUMENTS



LEAK TEST

Your Text

Unique features

- Extremely "COMPACT" construction
- Lasting vacuum pump
- Polycarbonate vacuum dessicator
- Stainless steel [GMP] CONSTRUCTION
- Vacuum level up to 22" of HG

Ideal instrument to check blister packs , strips for leak test. The old method of separate dessicator and vacuum pump is replaced with this new leak test apparatus in GMP Stainless steel body with phenomenal compactness the instrument is equipped with long lasting "INDIAN" vacuum pump which gives higher level of vacuum in shortest Possible time .Die pressed polycarbonate vacuum dessicator sustains vacuum of long time . A vacuum gauge is provided to indicate vacuum level which is connected to isolation valve to disconnect the source .thus avoiding the need of continuous running od vacuum pump . A 0 TO 5 minute timer helps the operator to do other work its amazing compactness saves lot of valuable place in in process laboratories.

Specifications :

Dessicator dia : 8" / 10" / 12"
Vacuum level : 22" Hg
Dessicator : polycarbonate

Rating : 5A 230 V single phase
Timer : 5 minute digital
Body : fully stainless steel



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinn's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com