



## MAGNETIC STIRRER (MODEL NO.BTI-33)



### SPECIFICATION

<b>Make</b>	<b>:</b>	<b>BTI</b>
<b>Top covers</b>	<b>:</b>	<b>S.S.304</b>
<b>Max Stirring Cap</b>	<b>:</b>	<b>2 MLH</b>
<b>Speed Range</b>	<b>:</b>	<b>100 to 1000 RPM</b>
<b>Temp. Range</b>	<b>:</b>	<b>40 °C to 120 °C</b>
<b>Power Supply</b>	<b>:</b>	<b>230 Volts, 50 H<sub>2</sub></b>

**Construction** : **Top cover made of S.S 304 and body made of M.S.**

**Powder coated.**

**Temperature controlled by : Energy Regulator**

