



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com



ORBITAL SHAKING INCUBATOR (HEATING)



Size	Capacity (Flasks)
a)	8x250ml.
b)	16x250ml.
c)	8x500ml.
d)	16x500ml.

ORBITAL SHAKING INCUBATOR (HEATING)

- Double Walled Inner S.S. & Outer Body M.S. Powder Coated.
- Temp.Range:5 °C above Room Temperature to max 60°C.
- Temp. Controlled by: Microprocessor PID Temp Controller
- RPM Range:50 to 200 RPM
- Variable Speed controlled by Digital Display with PMDC Motor complete with stainless steel tray Lotus shaped spring type clamp but without conical Flask.

SUNSHINE INSTRUMENTS