



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com



OIL BATH WITH DIGITAL CONTROLLER



SPECIFICATION

Construction: Double walled outer body made of mild steel powder coated / stainless steel & Inner body made of stainless steel 304

Special class glass wool insulation to prevent loss of heating

Special Teflonwiser oil heater made of S.S 304

Temperature controller by PID microprocessor controller

Accuracy +/- 2 °C

Temperature range 50 °C to 300 °C

Cap inLtr	Inner working size of Oil bath Size in inches	Rating watts
1 Ltr	8" x 6"	1000
2 Ltr	8"x8"	1200
3 Ltr	9" x 9"	1200
5 Ltr	10"x10"	1500
10 Ltr	12"x12"	2000

SUNSHINE INSTRUMENTS