



**IMPORTERS / EXPORTERS & SUPPLIERS :-**

# SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com



## ORBITAL SHAKER: (LAB USE)



### ORBITAL SHAKER: (LAB USE)

The Body made of M.S & painted with epoxy powder coated. The tray made of S.S. shaking speed: 20 to 200 orbital strokes per minute & the orbit diameter of 25mm supplied with S.S. clamp for conical flask.

Size	Size in MM	With analogue speed controller	With digital speed indicator
a) 8 x 250ml	325 X 325		
b) 16 x 250ml	425 X 425		
c) 8 x 500ml	325 X 325		
d) 16 x 500ml	425 X 425		
e) 4 x 1000ml	325 X 325		
f) 8 x 1000ml	425 X 425		

Optional Accessories:

Optional: Electronic Digital Timer ( 0 to 999 Minutes) for above

# SUNSHINE INSTRUMENTS