

## Refractometer PAL-1



The PAL-1 works perfectly in measuring almost any fruit juice, food, or drink, such as soup, sauce, ketchup, tomato sauce, low-sugar jam, or marmalade. PAL series comes with Data Transmission Function with NFC (Near Field Communication) technology. Data history can be accessed by bringing PAL-NFC to any Android devices or PC-linked USB NFC Reader.

### Specifications

Model	PAL-1
Cat.No.	3810
Range	Brix 0.0-53.0 %
Resolution	Brix 0.1%
Accuracy	±0.2 %



Shop



Quotation Request

### Accessories

- PAL-CASE : RE-39409
- JORDAN for PAL(Silicon Cover) : RE-39414
- STRAP : RE-39410
- SUPER EXTRACTOR : RE-29401
- 10% Sucrose (±0.03%) : RE-110010
- 20% Sucrose (±0.03%) : RE-110020
- 30% Sucrose (±0.03%) : RE-110030
- 40% Sucrose (±0.04%) : RE-110040
- 50% Sucrose (±0.05%) : RE-110050
- Battery Case Cover With O-Ring : RE-38163
- AC Battery Charger : RE-34802
- \* Calibration Certificate : Contact an ATAGO representative for details.