

## Pilot 100 Vacuum Controller

- ◆ One-year warranty
- ◆ Start and leave with evaporation detect and auto adjustments.
- ◆ GUI control panel, what you see is what you get.
- ◆ Leave your glove on, Pilot can still read your taps.
- ◆ Can be installed in every vacuuming sources.
- ◆ Suitable for laboratory or industrial environments.
- ◆ CE, UKCA, RoHS certificated



### Pilot 100 Vacuum Controller Specification

- ▶ **Model: Pilot 100**
- ▶ Power input: DC24V, 1A
- ▶ Resolution: 1 digit
- ▶ Pressure display: mbar, hPa, torr
- ▶ Display: 7" touch screen
- ▶ Measuring range: 1100 ~ 0 mbar abs.
- ▶ Measurement accuracy:  $\pm 2$  mbar (after adjustment, constant temp.)
- ▶ Setting range: 1013 ~ 0 mbar
- ▶ Housing material: Plastic
- ▶ Venting valve: Yes
- ▶ Auto-cleaning: Yes
- ▶ Power: 15 W
- ▶ Hose barb: ID 8 ~ ID10
- ▶ Safety: Over-pressure Protection, Error Message
- ▶ Dimension ( L x W x H ): 21 x 20 x 18 cm
- ▶ Net weight: 2.2 Kg