



Oil Free Compressor



Features

- ▶ **No pollution, maintenance free**
Rocker series pumps are driven by piston, without the need of lubricant, regular oil changes and maintenance; with no oil pollution.
- ▶ **Quiet and low vibration**
Driven directly by motor with no additional belt-driven transmission; the quality of vibration-proof assembly makes Rokcer series pumps run at the lowest noise level among all other equivalent pumps.
- ▶ **Thermal protection device**
Every motor of Rocker series pumps has a built-in thermal protection device to shut off the pump automatically when overheated and then resume working when the temperature cools down.
- ▶ **Air outlet filter device**
Every Rocker series compressor is equipped with filter cartridge at the air outlet to filter particle, moisture and supply clean air.
- ▶ **International certification**
CE certification
CSA certification
IP 30 protection
RoHS certification
- ▶ **Two-Year Warranty**

Applications

- ▶ Aeration
- ▶ Pressure filtration

Ordering Information

167320-11(22)

Rocker 320, Compressor
AC110V, 60Hz (AC220V, 50Hz)

167420-11(22)

Rocker 420, Compressor
AC110V, 60Hz (AC220V, 50Hz)



Rocker 320



Rocker 420

Specification

Model	Rocker 320		Rocker 420	
Voltage / Freq.	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz
Max. power	90W	85W	147W	125W
Max. current	0.9A	0.4A	1.5A	0.6A
Max. pressure (w/o moisture trap)	100 psi		100 psi	
Max. pressure (w. moisture trap)	80 psi		80 psi	
Max. flow rate	14 L/min	12 L/min	27 L/min	25 L/min
Moisture trap / Pressure regulator	Yes		Yes	
Motor rotation speed	1700 RPM	1400 RPM	1700 RPM	1400 RPM
Horse power	1/8 HP		1/6 HP	
Noise level	50 dB		52 dB	
Hose barb	ID8		ID8	
Safety	Thermal Protection		Thermal Protection	
Net weight	4.1 Kg		5.2 Kg	
Dimension (LxWxH)	27.2 x 13.6 x 16.5 cm		29.8 x 14.3 x 17.2 cm	