



**IMPORTERS / EXPORTERS & SUPPLIERS :-**

# SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com



## ROTARY SHAKER (Heavy Duty)



### ROTARY SHAKER (Heavy Duty)

Variable speed P.M.D.C. Geared motor continues duty . RPM

20 to 250 with S.S. tray lotus clamp type. Body made of Galivanted sheets powder coated.

Supplied with S.S clamp

Size & No of Flasks	With analogue sped controller	With digital speed indicator
8x 250ml		
16x 250ml		
25x 250ml		
36x 250ml		
48x 250ml		
64x 250ml		

Size & No of Flasks	With analogue sped controller	With digital speed indicator
8x 500ml		
16x 500ml		
25x 500ml		
36x 500ml		
48x 500ml		
64x 500ml		

Optional Accessories:

Optional: Electronic Digital Timer ( 0 to 999 Minutes) for above

# SUNSHINE INSTRUMENTS