



## **BRIX REFRACTOMETER SCM-1000**

The Smart Digital Refractometer that is always ready.

### **PRODUCT DESCRIPTION**

Range: Brix 0.0 ~ 55.0%

Resolution: Brix 0.1%

Accuracy: Brix  $\pm 0.2\%$

ATC Function included (Temperature Compensation 10 ~ 60°C )

Temperature Range: 0 ~ 70°C

Within 3 Seconds Measurement time

Auto shut-off function after 2 minutes

Automatic screen brightness control

Data Store: Up to 20

Usage Environment: 10 ~ 40 °C

Idle Sample Volume: 0.5ml

Unit: Water resistant (splash proof)

Power Source 1.5V x 2(AAA alkaline)

Size: 58(W) x 36(D) x 122(H)mm

weight: 149g (including batteries and detachable cover)

