



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com



SOXHLET EXTRACTION HEATING MANTLE TYPE (Supplied without glasspart & without clamp and boshead)

Fiber glass hand knitted heating surface having Nichrome Wire heating elements. Body made of M.S. sheet duly finished with powder coating with energy regulators. Useful for quick heating which save extraction time.

Size	Rating Watts	No. of Holes	No. of Regulators
250 ml	450	3	3
500 ml	600	3	3
1 Lt.	900	3	3
250 ml	600	4	4
500 ml	800	4	4
1 Lt.	1000	4	4
250 ml	900	6	6
500 ml	1200	6	6
1 Lt.	1800	6	6

(Note: Warranty Does Not Cover for Nett Coil)

SPARE PARTS FOR ABOVE

THREE FINGURE CLAMP - Rs.380/- EACH

BOSHEAD HEAVY DUTY – Rs.220/- EACH

