



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS



Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com

SEROLOGICAL WATER BATH

Double Walled filled with high grade glass wool insulation between the two chambers.
Outer Chamber made of M.S with powder coating and Inner chamber made of Stainless Steel.
Temp. Controller by thermostatic Temp. Range 5°C above ambient to 60°C.
Temp: Accuracy $\pm 2^{\circ}\text{C}$
(Supplied without test tube racks)

Inner working size of Water bath Size	Inner Dimension in Inches	Outer M.S, Inner S.S	Complete S.S.	Pyramidal Lid
2 Rack	10X7X6			
4 Rack	14X10X6			

Optional Accessories for Water bath:

a) Digital Temperature indicator cum controller



SUNSHINE INSTRUMENTS