



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS




Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com

Tablet Hardness Tester



Tablet Hardness Tester

Features

It consists of a Force gauge & a Slide mounted on a bracket Side is operated by handwheel
Force gauge can be checked for accuracy separately.

Traceable certificate provided if required

Zero adjust for spot calibration.

Initial holding force is accounted for unlike Nlonsanto.

Red Pointer retains the maximum reading'

Acrylic dome cover prevent broken table pieces flying of all over.