



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com



ULTRA VOILET INSPECTION CABINET (U.V.CABINET)

ULTRA VOILET INSPECTION CABINET (U.V.CABINET):

3 Tubes
Useful for viewing paper and thin layer chromatograms G.I.epoxy powder coated cabinet. With enclose long, short wave and white source.

UV CABINET

TECHNICAL SPECIFICATION

MOC:MILD STEEL, POWDER COATED.

SIZE: 3 TUBES

VOEW PIECE: ALUMINUM CASTING WITH EYE PROTECTION

LIGHT SOURCE: 3 NOS

A) SHORT WAVE-254NM-1NO.

B) LONG WAVE-365NM.

C) VISIBLE LIGHT

TUBE SIZE: 12"LONG.

WORKING: 3 INDIVIDUAL TOGGLE SWITCH PROVIDED FOR EACH LIGHT.

