



SIZES IN RECTANGULAR MODEL		
S.No.	SIZE IN CM.	No. of Shelves
1.	20 x 20 x 30 Cm.	2
2.	30 x 30 x 40 Cm.	2
3.	35 x 35 x 50 Cm.	3
4.	45 x 45 x 45 Cm.	3
5.	45 x 45 x 60 Cm.	3
6.	50 x 50 x 75 Cm.	4

Vacuum Oven (RECTANGULAR) SI 53

- Vacuum Ovens Are Double Walled.
- Inside Chamber Is Made Of Thick Stainless Steel 304 Quality.
- Out Side Chamber Made Of Mild Steel Sheet Pre Greased, Pretreated For Rust Proofing And Powder Coated.
- High Grade Glass Wool Insulation Between The Chambers To Avoid Heat Loss. Uniform Drying Under Vacuum.
- Temp. Ranges From 50°C To 200°C With $\pm 1^\circ\text{C}$ Accuracy.
- Vacuum Gauge of 760mm hg Displayed On Analog Guage.
- Temp. Controlled By Digital Temp. Indicator Cum Controller.
- Vacuum Chamber Is Wrapped With Set Of High Quality Heaters For Quick And Uniform Heating.
- Door Has A Thick Tempered Toughned Glass Window. Perfectly Sealed With Silcon Rubber Gasket For Observation Without Disturbing Thermal Condition.
- Proper Door Alignment To Ensure Excellent Vacuum Performance And Operators Safety.
- Positive Screw To Tighten The Lid To Obtain Perfect Seal.
- Fitted With Diaphragm Valves And Vacuum Guage.
- Control Panel With Temp. Controller, Switches And Indicating Lamps Provided On The Front Panel.
- Suitable For 230V, Ac, 50hz.
- Available In Rectangular Chamber Also.
- Vacuum Control And Relief Valves Located On Front Of Unit For Convenient Monitoring.



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com