



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com

57 WATER BATH Thermostatic- Rectangular

Double Walled filled with high grade glass wool insulation between the two chambers. Outer Chamber made of M.S with powder coating and Inner chamber made of Stainless Steel. Temp. Controlled by thermostatic. Temp. Range 5°C above ambient to 99°C. Supplied with 7.5 cms concentric Ring.

Inner working size of Water bath Size in inches	Rating watts	No. of Holes	Outer M.S.BODY Price Rs./Each	Complete S.S.BODY Price Rs./Each
35 x 25 x 10 cm	1000	6 Hole		
40 x 30 x 10 cm	1500	12 Hole		

58 WATER BATH-(Round) Thermostatic (Without Lid & Ring)

Double walled construction inner body S.S 304 outer body M.S Powder coated /complete S.S 304. Glass wool insulation between two walled. Temperature range 5°C above room temperature to up to 99°C. Temperature controlled by hydraulic thermostat . Accuracy +/- 3°C to 4°C. C [for better accuracy take water bath with optional accessories digital controller.] Supplied without ring & lid

Cap. In Lt.	Inner working size of Water bath Size in Inches	Rating watts	Outer M.S BODY & Inner S.S Price Rs./Each	Complete S.S.BODY Price Rs./Each
1 Ltr	6 x 8	1000		
2 Ltr	8 x 8	1000		
3 Ltr	9 x 9	1500		
5 Ltr	10 x 10	1500		
10 Ltr	12x12	2000		

Optional Accessories for Water bath:

a) Micro processor PID controller(for accuracy +/- 1 deg)

