



IMPORTERS / EXPORTERS & SUPPLIERS :-
SUNSHINE INSTRUMENTS



Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com

WATER BATH-(Round) Thermostatic (Without Lid & Ring)

Double walled construction inner body S.S 304 outer body M.S Powder coated /complete S.S 304. Glass wool insulation between two walled. Temperature range 5°C above room temperature to up to 99°C. Temperature controlled by hydraulic thermostat . Accuracy +/- 3°C to 4°C. C [for better accuracy take water bath with optional accessories digital controller.] Supplied without ring & lid

Cap. In Lt.	Inner working size of Water bath Size in Inches	Rating watts	Outer M.S BODY & Inner S.S	Complete S.S.BODY
1 Ltr	6 x 8	1000		
2 Ltr	8 x 8	1000		
3 Ltr	9 x 9	1500		
5 Ltr	10 x 10	1500		
10 Ltr	12x12	2000		

Optional Accessories for Water bath:

a) Micro processor PID controller(for accuracy +/- 1 deg)

