



OTOSCOPES

EYE, EAR, NOSE & THROAT



22820
PocketScope™ Otoscope
with "AA" Handle



21504
Insufflation Bulb

PocketScope™ Otoscope with Throat Illuminator

Small, convenient, and easy to use with bright, white halogen light and a built-in throat illuminator.

- Halogen light for true tissue color and consistent, long-lasting illumination
- Simply detach the otoscope head to use the throat illuminator or all-purpose penlight
- Fiber optics provide cool light with no reflections or obstructions
- Wide-angle viewing lens allows for instrumentation under magnification
- Sealed system for pneumatic otoscopy

22820 PocketScope Otoscope with "AA" Handle

22821 PocketScope Otoscope with "AA" Handle and Soft Case

22822 PocketScope Otoscope with "AA" Handle and Hard Case

22800 PocketScope Otoscope with Rechargeable Handle

22811 PocketScope Otoscope with Rechargeable Handle and Soft Case

22812 PocketScope Otoscope with Rechargeable Handle and hard Case

22801 Pocket Otoscope with Throat Illuminator, Convertible Handle, Batteries and Desk Charger with Soft Case

22802 Pocket Otoscope with Throat Illuminator, Convertible Handle, Batteries and Desk Charger with Hard Case

21111 PocketScope Otoscope (head only)

PocketScope™ Otoscope with Throat Illuminator: Accessories

21504 Insufflation Bulb with Tip for Diagnostic Otoscope

72600 2.5 V Rechargeable Battery for PocketScope

03400-U 2.5 V Halogen Lamp for PocketScope Instruments

211026-502 Otoscope Illuminator Assembly

72830 AA Handle with Batteries

72800 Convertible Handle with Rechargeable Battery