

TECHNICAL DATA SHEET

ATAGO REFRACTOMETER MASTER 53M

Metal Body Version

DESCRIPTION

The Atago Refractometer MASTER-53M is a good alternative choice when the water-resistant feature is not a priority, and the temperature compensation will be carried out manually. The Atago Refractometer MASTER-53M is made from metal and can be used anytime, anywhere. This digital refractometer measuring device has also been equipped with a sensitive sensor so that it is able to respond to measurement results in a short time.

FEATURES & BENIFITS

- Made from metal
- Simple operation – apply sample, press start, and reading is displayed in seconds
- No more reading errors – digital display for easy readings
- Sleek new design
- Automatic sample distribution



SPECIFICATIONS

The Atago Refractometer MASTER-53M is made from metal.

| | |
|----------------------|--------------------------|
| Model: | MASTER-53M |
| Cat. No: | 2353 |
| Range: | Brix: 0.0 to 53.0% |
| Minimum Scale: | 0.02% |
| Accuracy: | ±0.2% |
| Repeatability: | ± 0.1% |
| Dimensions & Weight: | 3.2 x 3.4 x 16.8cm, 130g |

ORDERING INFORMATION

| Product Code | Packaging |
|--------------|-----------|
| 53M Atago | Each |