

Digital Hand-held "Pocket" Refractometer PAL-3


[Enlargement](#)


Quantity [Request Quotation](#)

[about the icons](#)


Designed with a full Range (Brix 0.0 to 93.0%), the PAL-3 works perfectly for measuring almost any samples. This model has been modified to have a high measurement accuracy at $\pm 0.1\%$ Brix, the highest accuracy model in PAL series.

Specifications

Model	PAL-3
Cat.No.	3830
Range	Brix : 0.0 to 93.0 %
Accuracy	± 0.1 %
Temperature compensation	10 to 100 ° C

Options

- MAGIC™ (for PAL) : RE-39411
 - STRAP : RE-39410
 - PAL-CASE : RE-39409
 - PAL Silicone Cover : RE-39413
 - 10% Sucrose ($\pm 0.03\%$) : RE-110010
 - 20% Sucrose ($\pm 0.03\%$) : RE-110020
 - 30% Sucrose ($\pm 0.03\%$) : RE-110030
 - 40% Sucrose ($\pm 0.04\%$) : RE-110040
 - 50% Sucrose ($\pm 0.05\%$) : RE-110050
 - 60% Sucrose ($\pm 0.05\%$) : RE-110060
- * Calibration Certificate : Contact an ATAGO representative for details.

Patent for Invention Registration No. ZL200310103015.2 (China)
 Patent for Design Registration No. ZL03303431.1 (China), 089244 (Taiwan)
 Patent Granted in countries around the world.