

Wrist Action Shaker



Product Code : LSI-SM-WAS-01



(An ISO 13485: 2016 Certified Co.)

Wrist Action Shaking Machine

It has a heavy cast iron base with rugged construction and two side arms for holding flasks. Fitted with the heavy-duty motor with variable speed can hold flask from 100 to 1000ml. Capacity to Hold 2 Flask, 4 Flask or 8 flask.

Key Features

- ✓ Body: Made of Mild Steel Powder coated
- ✓ Motor: Heavy duty geared DC motor
- ✓ Speed : Variable Speed
- ✓ Shaking Angle: 0 to 15°
- ✓ Speed: 20 to 450 RPM
- ✓ Max Load: 10 lb.
- ✓ Clamps: 2 Clamps to hold four Vessels
- ✓ Oscillations: Max, 450 per minute
- ✓ Timer: 1 to 60 Minutes Timer or Continuous running
- ✓ Unit Dimensions: 17.5" L x 17.5" W x 12" H
- ✓ Unit Weight: @ 15 Kgs.
- ✓ Mounted on: Heavy duty Vacuum Cups
- ✓ Power Supply: 220 V Single Phase

Lab Solution India

Mfr. & Suppliers of Laboratory Equipments
(CE, ISO 13485: 2016 Certified Company)

Office: Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

Whats app :- +9182957-32110; Mobile : +917056732110; Landline No.: 0171-2690051

Email: labsolution.india@gmail.com; Web Address: www.laboratoryinstrument.co.in , www.ishmayra.com

"We aim at Customer Delight rather than just Customer Satisfaction"