



SUNSHINE INSTRUMENTS

83, CHINNY'S CHAMBER, DR, NANJAPPA ROAD,
COIMBATORE-641009 (T.N) INDIA.

Email: instrumentsindia2012@gmail.com

Web :- www.lutrontestinginstruments.com

Ph - +91 9360023357 / 0422 2393357

"JAYANT" Test Sieves Wire Mesh Series Regular Used Specification 8" DIAMETER "A" CLASS WITH CC

B.S.S (410/1969) MESH NO.	A.S.T.M (11-70) MESH NO.	Indian Std (469/1972)	Microns
4	5	4.000 mm	4000
5	6	3.350 mm	3353
6	7	2.800 mm	2812
7	8	2.360 mm	2411
8	10	2.000 mm	2057
10	12	1.700 mm	1680
12	14	1.400 mm	1405
14	16	1.180 mm	1204
16	18	1.000 mm	1003
18	20	0.850 mm	850
22	25	0.710 mm	710
25	30	0.600 mm	600
30	35	0.500 mm	500
36	40	0.425 mm	420
44	45	0.355 mm	355
52	50	0.300 mm	300
60	60	0.250 mm	250
72	70	0.212 mm	210
85	80	0.180 mm	180
100	100	0.150 mm	150
120	120	0.125 mm	120
150	140	0.106 mm	105
170	170	0.090 mm	90
200	200	0.075 mm	75
240	230	0.063 mm	63
300	270	0.053 mm	53
350	325	0.045 mm	45
400	400	0.037 mm	37
500	500	0.025 mm	25
LID	AND RECEIVER	SET	

Stockist For :- FLUKE / RAYTEK / LUTRON / HTC / MECO / METRAVI / ROCKER / HM / EUTECH / LOVIBOND /
SHIMADZU / WHATMAN / REMI / HAMILTON / ERMA / ATAGO / MILWAUKEE / ZEAL / OPTICGLASS / LEIMCO
GSTN NO. – 33ADQPV5864Q1ZC IE Code – 3212008300 PAN:ADQPV5864Q

THIS IS A COMPUTER GENERATED COPY AND DOSE'NT REQUIRE ANY STAMP OR SIGN

Prices are subject to change without prior notice & based on Market Fluctuation. Page 1

Subject To Coimbatore Jurisdiction.

THE INSTRUMENTS FOR FUTURE