



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com



CAT NO J - 95 GYRATORY

GYRATORY SIEVE SHAKERS : "JAYANT" BRAND GYRATORY is heavy duty floor model. The Shaker is suitable for heavy materials; powders or Granules.

Specifications :


- | Fix Gyrotory Principle, Circular & Tapping from bottom of rotating plate.
- | All moving parts, gear's works in machine oil No. 40.
- | Standard Model suitable for Jayant Test Sieves 8" x 2" HT & with Adapter can be suitable for 12" x 4" HT Sieves.
- | Crompton 1/4 HP Motor 230 V, 50 Hz 180 Watts.
- | Timer Box, Separate unit built in with 3 pin socket, On/Off switch Pilot Lamp & 3 Core Cable with 3 pin plug 15 Amps & Mechanical Timer 10 to 60 minutes, minimum 10 minutes to set required time for sieve analysis.

Capacity : Maximum 8" x 2" HT 7 Sieves Plus U2

Maximum 12" HT x 4" HT 4 Sieves Plus U2

Weight : Unit : 40 kg. Shipping : 55 kg

Note : Rates without Test Sieves & Lid / Receiver Wooden case with Nut / Bolt packing.



SUNSHINE INSTRUMENTS