



IMPORTERS / EXPORTERS & SUPPLIERS :-

SUNSHINE INSTRUMENTS

Deals in :- Industrial / R&D / Scientific / Analytical / Diagnostics / Accessories

83, Chinny's Chamber, Dr. Nanjappa Road, Coimbatore - 641 018. (T.N)

Ph : +91 422 2393357, +91 93600 23357 E-mail : instrumentsindia2012@gmail.com

Website : www.lutrontestinginstruments.com



The Comparator is meant for visual determination of the colour of lubricating oils, heating oils, diesel fuel oils and petroleum waxes.

The instrument has been designed to conform specifications - US ASTM D 1500-64 and IS 1448:12

The standard colours are available in two discs

Disc A : 0.5, 1.0, 1.5, 2.0, 2.5, 3.0, 3.5, 4.0 (Price list code : 13.1.1)

Disc B : 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0 (Price list code : 13.1.2)

Procedure: The sample under test is taken in a 33mm dia test tube and is placed in the comparator. The disc is rotated until a match is obtained. Direct reading is then made in graded scale.

Size :12.5 cm x 12.5 cm x 10.5 cm (with wooden carry case)

Weight : 525 gms (with wooden carry case)